



## FREESTANDING BLOODBANK GLASS DOOR (360x250ml bags)



- ✓ **Innovative Controller with LCD Display**  
Fully meets MHRA guidelines
- ✓ **Heart Temperature Monitoring Software**  
Records air and load temperatures constantly, holding a full 7 day alarm history
- ✓ **+2°C to +6°C Load Temperature (bag temp) Range**
- ✓ **USB Slot Control Panel**  
A USB stick can be plugged into control panel, allowing user to download temperature information as often as required
- ✓ **Internal LED Light**
- ✓ **Lockable Castors**  
for ease of movement
- ✓ **Self Closing Door**  
Doors will close when positioned anything less than 90° open, and self lock for maximum security
- ✓ **Auto Defrost Feature**
- ✓ **Fully Removable, Vented Drawers**  
Easy-to-use drawer system offers maximum durability and allows cool air to follow correct path inside unit.
- ✓ **Easy Clean Interior**

SKU CODE:

PRODUCT CODE:

DESCRIPTION:

444440362

BBK360G

Freestanding Bloodbank Glass Door (360x250ml bags)



## FULL SPECIFICATION – BBK360G FREESTANDING BLOODBANK GLASS DOOR (360x250ml bags)

<b>Description</b>	Freestanding
<b>SKU Code</b>	444440362
<b>GDPA Code</b>	BBK360G
<b>No. of blood bags (250ml fill)</b>	360
<b>No. of blood bags (450ml fill)</b>	180
<b>Capacity (litres)</b>	350
<b>Working Temperature</b>	+4°C
<b>Load Temperature (bag temp)</b>	+2°C to +6°C
<b>Door</b>	Single
<b>Door Type</b>	Glass
<b>Door Construction</b>	Double Glazed
<b>Self Closing Door</b>	Yes
<b>Unit External Width (mm)</b>	680
<b>Unit External Depth (mm)</b>	910 (includes handle depth of 45mm)
<b>Unit External Height (mm)</b>	1935
<b>Product Weight (kg)</b>	216
<b>Packed Width (mm)</b>	740
<b>Packed Depth (mm)</b>	875
<b>Packed Height (mm)</b>	2150
<b>Packed Weight (kg)</b>	228
<b>Internal Dimensions (mm) (WxDxH)</b>	550 x 640 x 1000
<b>No. of Drawers</b>	5
<b>Lock Type</b>	Digital (with 2 key fobs)
<b>Digital display (LCD)</b>	Yes
<b>Internal Light (LED)</b>	Yes
<b>Core Temperature (audible/digital)</b>	Yes
<b>Air Temperature Alarm</b>	Yes
<b>Door Open Alarm</b>	Yes
<b>Mains Failure Alarm</b>	Yes
<b>Remote Alarm contacts</b>	Yes
<b>Fan Fail Alarm</b>	Yes
<b>Back-up Battery Alarm (visible)</b>	Yes
<b>Compressor Battery Back-up</b>	Yes
<b>Load, Air &amp; Defrost Sensor Fail Alarm</b>	Yes
<b>Controller Battery Back-up</b>	Yes
<b>USB</b>	Yes
<b>Chart Recorder</b>	Yes
<b>Ambient working Temperature</b>	Max: +43°C, Min: +10°C
<b>Castors</b>	Yes (4, 2 lockable at front)
<b>Adjustable Feet</b>	-
<b>Fixed Feet</b>	-
<b>Fan Type</b>	Axial
<b>Automatic Defrost</b>	Yes
<b>Dual Refrigeration</b>	No
<b>Refrigerant</b>	R600a
<b>Medical Devices Certification</b>	Yes
<b>BS4376 Compliant</b>	Yes
<b>Voltage</b>	220v-240v 50Hz/60Hz