



FREESTANDING BLOODBANK SOLID DOOR (864x250ml bags)



- ✓ **Innovative Controller with LCD Display**
Fully meets MHRA guidelines
- ✓ **Heart Temperature Monitoring Software**
Records air and load temperatures constantly, holding a full 7 day alarm history
- ✓ **+2°C to +6°C Load Temperature (bag temp) Range**
- ✓ **USB Slot Control Panel**
A USB stick can be plugged into control panel, allowing user to download temperature information as often as required
- ✓ **Internal LED Light**
- ✓ **Lockable Castors**
for ease of movement
- ✓ **Self Closing Door**
Doors will close when positioned anything less than 90° open, and self lock for maximum security
- ✓ **Auto Defrost Feature**
- ✓ **Fully Removable, Vented Drawers**
Easy-to-use drawer system offers maximum durability and allows cool air to follow correct path inside unit
- ✓ **Easy Clean Interior**

SKU CODE:

PRODUCT CODE:

DESCRIPTION:

444440884

BBK864S

Freestanding Bloodbank Solid Door (864x250ml bags)



FULL SPECIFICATION – BBK864S FREESTANDING BLOODBANK SOLID DOOR (864x250ml bags)

| | |
|---------------------------------------|------------------------------|
| Description | Freestanding |
| SKU Code | 444440884 |
| GDPA Code | BBK864S |
| No. of blood bags (250ml fill) | 864 |
| No. of blood bags (450ml fill) | 432 |
| Capacity (litres) | 960 |
| Working Temperature | +4°C |
| Load Temperature (bag temp) | +2°C to +6°C |
| Door | Double |
| Door Type | Solid |
| Self Closing Door | Yes |
| Door Construction | Zintec Steel |
| Unit External Width (mm) | 1340 |
| Unit External Depth (mm) | 850 |
| Unit External Height (mm) | 1980 |
| Product Weight (kg) | 370 |
| Packed Width (mm) | - |
| Packed Depth (mm) | - |
| Packed Height (mm) | - |
| Packed Weight (kg) | 400 |
| Internal Dimensions (mm) (WxDxH) | 1225 x 640 x 1225 |
| No. of Drawers | 12 |
| Lock Type | Digital (with 2 key fobs) |
| Digital display (LCD) | Yes |
| Internal Light (LED) | Yes |
| Core Temperature (audible/digital) | Yes |
| Air Temperature Alarm | Yes |
| Door Open Alarm | Yes |
| Mains Failure Alarm | Yes |
| Remote Alarm contacts | Yes |
| Fan Fail Alarm | Yes |
| Back-up Battery Alarm (visible) | Yes |
| Compressor Battery Back-up | Yes |
| Load, Air & Defrost Sensor Fail Alarm | Yes |
| Controller Battery Back-up | Yes |
| USB | Yes |
| Chart Recorder | Yes |
| Ambient working Temperature | - |
| Castors | Yes (4, 2 lockable at front) |
| Adjustable Feet | - |
| Fixed Feet | - |
| Fan Type | Axial |
| Automatic Defrost | Yes |
| Dual Refrigeration | Yes |
| Refrigerant | R134a |
| Medical Devices Certification | Yes |
| BS4376 Compliant | Yes |
| Voltage | 220v-240v 50Hz/60Hz |