



UNDER-COUNTER LABORATORY FREEZER SOLID DOOR 102L WHITE



- ✓ **-18°C to -20°C Operating Temperature Range**
Conforms to MHRA guidelines
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Wire Shelves and Basket**
Enables optimum air flow inside unit

SKU CODE:

PRODUCT CODE:

DESCRIPTION:

444441803

ISU37C

Under-counter Laboratory Freezer Solid Door 102L White



FULL SPECIFICATION – ISU37C UNDER-COUNTER LABORATORY FREEZER SOLID DOOR 102L WHITE

Description	Under-counter
SKU Code	444441803
EAN Code	5034648493271
GDPA Code	ISU37C
Gross Capacity (L)	102
Net Capacity (L)	87
Unit External Width (mm)	495
Unit External Depth (mm)	620
Unit External Height (mm)	830
Product Weight (kg)	33
Packed Width (mm)	560
Packed Depth (mm)	655
Packed Height (mm)	850
Packed Weight (kg)	36.5
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	414 x 415 x 683
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	-18°C to -20°C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1050
Replacement Door Seal	No
Open Door Sensor	No
Lock	Yes
Adjustable Feet	1, front
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	2
Number of Shelf Positions	Fixed
Shelf Size	-
Shelf Material	Wire
Moulded Drawers	0
Adjustable Door Bins	0
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Manual
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R134a
CFC Free	Yes
HFC Free	No
HCFC Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	No
Insulation Thickness (mm)	35
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes