



FREESTANDING LABORATORY FRIDGE SOLID DOOR 475L WHITE



- ✓ **+2°C to +10°C Operating Range**
Conforms to MHRA guidelines
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Castors**
For ease of movement
- ✓ **Wire Shelves**
Enables optimum air flow inside unit

SKU CODE: PRODUCT CODE:

444441801 LR1607C

DESCRIPTION:

Freestanding Laboratory Fridge Solid Door 475L White



FULL SPECIFICATION – LR1607C FREESTANDING LABORATORY FRIDGE SOLID DOOR 475L WHITE

Description	Freestanding
SKU Code	444441801
EAN Code	5034648485696
GDPA Code	LR1607C
Gross Capacity (L)	475
Net Capacity (L)	441
Unit External Width (mm)	695
Unit External Depth (mm)	910
Unit External Height (mm)	1688
Product Weight (kg)	95
Packed Width (mm)	770
Packed Depth (mm)	950
Packed Height (mm)	1810
Packed Weight (kg)	116
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	550 x 670 x 1170
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	+2°C to +10°C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, R/H only
Door Open Depth (mm)	1596
Replacement Door Seal	No
Open Door Sensor	No
Lock	Yes
Adjustable Feet	Castors
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	5
Number of Shelf Positions	30
Shelf Size	5: 536 x 609 x 6
Shelf Material	Wire
Moulded Drawers	0
Adjustable Door Bins	0
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Automatic
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R134a
CFC Free	Yes
HFC Free	No
HCFC Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	No
Insulation Thickness (mm)	65
Insulation Type	Cyclopentane
Anti Microbial Handle	No