




COUNTERTOP LABORATORY FRIDGE SOLID DOOR 82L WHITE



- ✓ **+2°C to +10°C Operating Temperature Range**
Conforms to MHRA guidelines
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Wire Shelves and Basket**
Enables optimum air flow inside unit
- ✓ **Audible Warning Alarms**
Alerts user should internal temperature deviate from set parameters

SKU CODE: PRODUCT CODE:

444441799 LR207C

DESCRIPTION:

Countertop Laboratory Fridge Solid Door 82L White



FULL SPECIFICATION – LR207C COUNTERTOP LABORATORY FRIDGE SOLID DOOR 82L WHITE

Description	Countertop
SKU Code	444441799
EAN Code	5034648485689
GDPA Code	LR207C
Gross Capacity (L)	82
Net Capacity (L)	67
Unit External Width (mm)	502
Unit External Depth (mm)	540
Unit External Height (mm)	660
Product Weight (kg)	27.5
Packed Width (mm)	520
Packed Depth (mm)	595
Packed Height (mm)	675
Packed Weight (kg)	30.5
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	415 x 340 x 540
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	+2°C to +10°C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	976
Replacement Door Seal	No
Open Door Sensor	No
Lock	Yes
Adjustable Feet	No
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	2
Number of Shelf Positions	2
Shelf Size	2: 438 x 274 x 6
Shelf Material	Wire
Moulded Drawers	0
Adjustable Door Bins	0
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Automatic
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R134a
CFC Free	Yes
HFC Free	No
HCFC Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	No
Insulation Thickness (mm)	30
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes