




## UNDER-COUNTER LABORATORY FRIDGE SOLID DOOR 107L WHITE



- ✓ **+2°C to +10°C Operating Temperature Range**  
Conforms to MHRA guidelines
- ✓ **NEW Unique Anti Microbial Handle**  
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**  
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Wire Shelves and Basket**  
Enables optimum air flow inside unit
- ✓ **Audible Warning Alarms**  
Alerts user should internal temperature deviate from set parameters

SKU CODE:    PRODUCT CODE:

444441800    LR307C

DESCRIPTION:

Under-counter Laboratory Fridge Solid Door 107L White



## FULL SPECIFICATION – LR307C UNDER-COUNTER LABORATORY FRIDGE SOLID DOOR 107L WHITE

<b>Description</b>	Under-counter
<b>SKU Code</b>	444441800
<b>EAN Code</b>	5034648493264
<b>GDPA Code</b>	LR307C
<b>Gross Capacity (L)</b>	107
<b>Net Capacity (L)</b>	103
<b>Unit External Width (mm)</b>	495
<b>Unit External Depth (mm)</b>	620
<b>Unit External Height (mm)</b>	830
<b>Product Weight (kg)</b>	31.5
<b>Packed Width (mm)</b>	560
<b>Packed Depth (mm)</b>	655
<b>Packed Height (mm)</b>	850
<b>Packed Weight (kg)</b>	35
<b>Country of Origin</b>	UK
<b>Internal Dimensions (mm) (WxDxH)</b>	403 x 400 x 540
<b>External Finish</b>	White
<b>Internal Finish</b>	White
<b>Type of Plastic</b>	SRL800 (HIPS/ABS Mix)
<b>Temperature Range</b>	+2°C to +10°C
<b>Colour</b>	White/Green
<b>Door Type</b>	Solid
<b>Door Style</b>	Hinged, Reversible
<b>Door Open Depth (mm)</b>	1050
<b>Replacement Door Seal</b>	No
<b>Open Door Sensor</b>	No
<b>Lock</b>	Yes
<b>Adjustable Feet</b>	2, front
<b>Temperature Display on Panel</b>	Yes
<b>Temperature Alarm</b>	Yes
<b>Number of Shelves</b>	2
<b>Number of Shelf Positions</b>	4
<b>Shelf Size</b>	2: 427 x 374 x 6
<b>Shelf Material</b>	Wire
<b>Moulded Drawers</b>	0
<b>Adjustable Door Bins</b>	0
<b>Interior Light</b>	No
<b>Fan</b>	No
<b>Spark Free</b>	Yes
<b>Auto or Manual defrost</b>	Automatic
<b>Controller type</b>	Digital, External
<b>Climate Class</b>	N
<b>Refrigerant Type</b>	R134a
<b>CFC Free</b>	Yes
<b>HFC Free</b>	No
<b>HCFC Free</b>	Yes
<b>UK 3 Pin Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Voltage Supply</b>	230
<b>Aux Output</b>	No
<b>Insulation Thickness (mm)</b>	35
<b>Insulation Type</b>	Cyclopentane
<b>Anti Microbial Handle</b>	Yes