

Spark-Free Fridge-Freezer Range

Offering the flexibility of a fridge and freezer in the same unit, and maximising storage space, Lec Medical provides a comprehensive range of combination fridge-freezers in different sizes to suit your requirements.

Key Features

- ✓ Spark-free to avoid potential combustion
- ✓ External digital controller
- ✓ Lock and two keys for improved safety
- ✓ Antimicrobial handle

LabCool.de • Division of National Lab GmbH
Grambeker Weg 157 • 23879 Moelln • Germany
e-Mail: LabSales@NationalLab.com

Freestanding Laboratory
Fridge-Freezer Combi Solid
Door 263L



FRIDGE
2 to 10°C


FREEZER
-12 to -20°C



Specification

Freestanding Laboratory Fridge-Freezer Combi Solid Door 263L	
GDPA code	LSC263UK
SKU code	444440973
Gross capacity	263 (Fridge: 159, Freezer: 104)
Number of shelves	3 (Fridge)
Number of shelf positions	6
Refrigerant type	R600a
Temperature range	Fridge: 2°C to 10°C, Freezer: -12°C to -20°C
Voltage supply	230

Dimensions



Also Available...



Under-Counter Laboratory Fridge-Freezer Combi Solid Door 119L

LSC119UK / 444440972

Freestanding Laboratory Fridge-Freezer Combi Solid Door 324L

LSC324UK / 444441790



Fridge-freezers offer the optimal storage flexibility with two distinct temperature zones.

These units are a popular choice when space is limited and we offer a variety of sizes, from countertop units through to large freestanding alternatives.