



## FREESTANDING LABORATORY FREEZER SOLID DOOR 232L WHITE



- ✓ **-18°C to -25°C Operating Temperature Range**  
Desired temperature can be set by the user
- ✓ **NEW Advanced Controller**  
Records temperatures accurate to 0.1°C
- ✓ **NEW Unique Anti Microbial Handle**  
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**  
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**  
For ease of movement
- ✓ **Plastic Drawers and Flaps**  
Maximises space and offers protection to contents
- ✓ **Audible Warning Alarms**  
Alerts user should internal temperature deviate from set parameters

SKU CODE:    PRODUCT CODE:

444440970    LSF232UK

DESCRIPTION:

Freestanding Laboratory Freezer Solid Door 232L White



## FULL SPECIFICATION – LSF232UK FREESTANDING LABORATORY FREEZER SOLID DOOR 232L WHITE

|   |                                      |
|---|--------------------------------------|
| <b>Description</b>                      | Freestanding Laboratory Freezer      |
| <b>SKU Code</b>                         | 444440970                            |
| <b>EAN Code</b>                         | 5052263009707                        |
| <b>GDPA Code</b>                        | LSF232UK                             |
| <b>Gross Capacity (L)</b>               | 232                                  |
| <b>Net Capacity (L)</b>                 | 212                                  |
| <b>Unit External Width (mm)</b>         | 595                                  |
| <b>Unit External Depth (mm)</b>         | 634                                  |
| <b>Unit External Height (mm)</b>        | 1570                                 |
| <b>Product Weight (kg)</b>              | 60                                   |
| <b>Packed Width (mm)</b>                | 657                                  |
| <b>Packed Depth (mm)</b>                | 734                                  |
| <b>Packed Height (mm)</b>               | 1585                                 |
| <b>Packed Weight (kg)</b>               | 65                                   |
| <b>Country of Origin</b>                | UK                                   |
| <b>Internal Dimensions (mm) (WxDxH)</b> | 432 x 430 x 1300 (step: 190d x 255h) |
| <b>External Finish</b>                  | White                                |
| <b>Internal Finish</b>                  | White                                |
| <b>Type of Plastic</b>                  | SRL800 (HIPS/ABS Mix)                |
| <b>Temperature Range</b>                | -18°C to -25°C                       |
| <b>Colour</b>                           | White/Green                          |
| <b>Door Type</b>                        | Solid                                |
| <b>Door Style</b>                       | Hinged, Reversible                   |
| <b>Door Open Depth (mm)</b>             | 1185                                 |
| <b>Replacement Door Seal</b>            | No                                   |
| <b>Open Door Sensor</b>                 | Yes                                  |
| <b>Lock</b>                             | Yes                                  |
| <b>Adjustable Feet</b>                  | 2 Adjustable Feet, 2 Rear Rollers    |
| <b>Temperature Display on Panel</b>     | Yes                                  |
| <b>Temperature Alarm</b>                | Yes                                  |
| <b>Number of Shelves</b>                | 4 Drawers, 2 Fixed Flaps             |
| <b>Number of Shelf Positions</b>        | Fixed                                |
| <b>Shelf Size</b>                       | -                                    |
| <b>Shelf Material</b>                   | Wire Baskets                         |
| <b>Moulded Drawers</b>                  | 4 Solid                              |
| <b>Adjustable Door Bins</b>             | 0                                    |
| <b>Interior Light</b>                   | No                                   |
| <b>Fan</b>                              | No                                   |
| <b>Spark Free</b>                       | Yes                                  |
| <b>Auto or Manual defrost</b>           | Manual                               |
| <b>Controller type</b>                  | Digital, External                    |
| <b>Climate Class</b>                    | SN                                   |
| <b>Refrigerant Type</b>                 | R600a                                |
| <b>CFC Free</b>                         | Yes                                  |
| <b>HFC Free</b>                         | Yes                                  |
| <b>HCFC Free</b>                        | Yes                                  |
| <b>UK 3 Pin Plug</b>                    | Yes                                  |
| <b>Other Plug Types Available</b>       | Yes                                  |
| <b>Voltage Supply</b>                   | 230                                  |
| <b>Aux Output</b>                       | Yes                                  |
| <b>Insulation Thickness (mm)</b>        | 60                                   |
| <b>Insulation Type</b>                  | Cyclopentane                         |
| <b>Anti Microbial Handle</b>            | Yes                                  |