



UNDER-COUNTER PHARMACY FRIDGE SOLID DOOR 107L WHITE



- ✓ **+2°C to +8°C Operating Temperature Range**
Conforms to MHRA guidelines
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Wire Shelves and Basket**
Enables optimum air flow inside the unit

SKU CODE: PRODUCT CODE:

444441794 PE307C

DESCRIPTION:

Under-counter Pharmacy Fridge Solid Door 107L White



FULL SPECIFICATION – PE307C UNDER-COUNTER PHARMACY FRIDGE SOLID DOOR 107L WHITE

Description	Under-counter
SKU Code	444441794
EAN Code	5052263017948
GDPA Code	PE307C
Gross Capacity (litres)	107
Net Capacity (litres)	103
Unit External Width (mm)	495
Unit External Depth (mm)	620
Unit External Height (mm)	830
Product Weight (kg)	31.5
Packed Width (mm)	560
Packed Depth (mm)	655
Packed Height (mm)	850
Packed Weight (kg)	35
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	403 x 400 x 700
Temperature Range	+2°C to +8°C
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1050
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet (number & position)	1 adjustable on front
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	2
Number of Shelf Positions	4
Shelf Size (mm) (WxDxH)	2: 427 x 374 x 6
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Interior Light	No
Fan	Yes
Spark Free	No
Auto or Manual Defrost	Automatic
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R134a
CFC Free	Yes
HFC Free	No
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	35
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	-