




## FREESTANDING LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 263L WHITE



- ✓ **Operating Temperature Range of +2°C to +10°C for Fridge and -18°C to -25°C for Freezer**  
Desired temperature can be set by the user
- ✓ **NEW Advanced Controller**  
Records temperatures accurate to 0.1°C
- ✓ **NEW Unique Anti Microbial Handle**  
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **2 Fitted Locks as Standard (one on fridge, one on freezer)**  
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**  
For ease of movement
- ✓ **Wire Shelves**  
Enables optimum air flow inside unit
- ✓ **Plastic Baskets and Door Organisers**  
Maximises space and offers protection to contents
- ✓ **Audible Warning Alarms**  
Alerts user should internal temperature deviate from set parameters

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444440973	LSC263UK	Freestanding Laboratory Fridge-freezer Combi Solid Door 263L White



## FULL SPECIFICATION – LSC263UK FREESTANDING LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 263L WHITE

<b>Description</b>	Freestanding Laboratory Fridge-freezer Combi
<b>SKU Code</b>	444440973
<b>EAN Code</b>	5052263009738
<b>GDPA Code</b>	LSC263UK
<b>Gross Capacity (L)</b>	263 (Fridge: 159, Freezer: 104)
<b>Net Capacity (L)</b>	245
<b>Unit External Width (mm)</b>	540
<b>Unit External Depth (mm)</b>	600
<b>Unit External Height (mm)</b>	1680
<b>Product Weight (kg)</b>	58
<b>Packed Width (mm)</b>	597
<b>Packed Depth (mm)</b>	684
<b>Packed Height (mm)</b>	1700
<b>Packed Weight (kg)</b>	64
<b>Country of Origin</b>	UK
<b>Internal Dimensions (mm) (WxDxH)</b>	Fridge: 457 x 450 x 725, Freezer: 408 x 415 x 670 (step: 170d x 220h)
<b>External Finish</b>	White
<b>Internal Finish</b>	White
<b>Type of Plastic</b>	SRL800 (HIPS/ABS Mix)
<b>Temperature Range</b>	Fridge: +2°C to +10°C, Freezer: -18°C to -25 °C
<b>Colour</b>	White/Green
<b>Door Type</b>	Solid
<b>Door Style</b>	Hinged, Reversible
<b>Door Open Depth (mm)</b>	1110
<b>Replacement Door Seal</b>	No
<b>Open Door Sensor</b>	Yes
<b>Lock</b>	2 fitted locks (one on fridge, one on freezer)
<b>Adjustable Feet</b>	2 adjustable front feet, 2 fixed rear
<b>Temperature Display on Panel</b>	Yes
<b>Temperature Alarm</b>	Yes
<b>Number of Shelves</b>	3 (Fridge)
<b>Number of Shelf Positions</b>	6
<b>Shelf Size</b>	3: 308 x 448 x 4.5
<b>Shelf Material</b>	Wire
<b>Moulded Drawers</b>	2 Solid half (Fridge)
<b>Adjustable Door Bins</b>	2 full + 2 half
<b>Interior Light</b>	No
<b>Fan</b>	No
<b>Spark Free</b>	Yes
<b>Auto or Manual defrost</b>	Fridge: Automatic, Freezer: Manual
<b>Controller type</b>	Digital, External
<b>Climate Class</b>	SN
<b>Refrigerant Type</b>	R600a
<b>CFC Free</b>	Yes
<b>HFC Free</b>	Yes
<b>HCFC Free</b>	Yes
<b>UK 3 Pin Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Voltage Supply</b>	230
<b>Aux Output</b>	Yes
<b>Insulation Thickness (mm)</b>	60
<b>Insulation Type</b>	Cyclopentane
<b>Anti Microbial Handle</b>	Yes