




FREESTANDING FRIDGE-FREEZER COMBI SOLID DOOR 324L WHITE



- ✓ **Operating Temperature Range of +2°C to +10°C for Fridge and -18°C to -25°C for Freezer**
Conforms to MHRA guidelines
- ✓ **NEW Advanced Controller**
Records temperatures accurate to 0.1°C
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **2 Fitted Locks as Standard (one on fridge, one on freezer)**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**
For ease of movement
- ✓ **Wire Shelves**
Enables optimum air flow inside unit
- ✓ **Plastic Baskets and Door Organisers**
Maximises space and offers protection to contents

SKU CODE: PRODUCT CODE:

444441790 LSC324UK

DESCRIPTION:

Freestanding Fridge-freezer Combi Solid Door 324L White



FULL SPECIFICATION – LSC324UK FREESTANDING LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 324L WHITE

Description	Freestanding
SKU Code	444441790
EAN Code	5052263017900
GDPA Code	LSC324UK
Gross Capacity (L)	324 (Fridge: 190, Freezer: 134)
Net Capacity (L)	302
Unit External Width (mm)	600
Unit External Depth (mm)	600
Unit External Height (mm)	1870
Product Weight (kg)	70
Packed Width (mm)	657
Packed Depth (mm)	734
Packed Height (mm)	1885
Packed Weight (kg)	78
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	Fridge: 510 x 430 x 805, Freezer: 463 x 445 x 780 (step:250d x 260h)
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	Fridge: +2°C to +10°C, Freezer: -18°C to -25 °C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1190
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet	2 adjustable feet, 2 rear rollers
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	3 (Fridge)
Number of Shelf Positions	6
Shelf Size	3: 320 x 505 x 4.5
Shelf Material	Wire
Moulded Drawers	2 Solid half (Fridge)
Adjustable Door Bins	2 full + 4 half
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Fridge: Automatic, Freezer: Manual
Controller type	Digital, External
Climate Class	SN
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	60
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes