



UNDER-COUNTER PHARMACY FRIDGE SOLID DOOR 151L WHITE



- 
+2°C to +8°C Operating Temperature Range
 Conforms to MHRA guidelines
- 
NEW Advanced Controller
 Records temperatures accurate to 0.1°C
- 
NEW Unique Anti Microbial Handle
 Helps prevent the growth of bacteria and minimise the potential for cross contamination
- 
Fitted Lock as Standard
 Maintains security and ensures the unit can only be accessed by approved personnel
- 
Fitted Rollers
 For ease of movement
- 
Wire Shelves and Basket
 Enables optimum air flow inside the unit
- 
Audible Warning Alarms
 Alerts user should internal temperature deviate from set parameters

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444440940	PSR151UK	Under-counter Pharmacy Fridge Solid Door 151L White



FULL SPECIFICATION – PSR151UK UNDER-COUNTER PHARMACY FRIDGE SOLID DOOR 151L WHITE

Description	Under-counter
SKU Code	444440940
EAN Code	5052263009400
GDPA Code	PSR151UK
Gross Capacity (litres)	151
Net Capacity (litres)	141
Unit External Width (mm)	595
Unit External Depth (mm)	610 (w/o handle), 635 (with handle)
Unit External Height (mm)	845
Product Weight (kg)	37
Packed Width (mm)	657
Packed Depth (mm)	694
Packed Height (mm)	872
Packed Weight (kg)	41
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	500 x 450 x 705 (step: 165d x 185h)
Temperature Range	+2°C to +8°C
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1165
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	2 Adjustable Front Feet, 2 Rear Rollers
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	3
Number of Shelf Positions	6
Shelf Size (mm) (WxDxH)	3: 400 x 520 x 6
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Interior Light	No
Fan	Yes
Spark Free	No
Auto or Manual Defrost	Automatic
Controller type	Digital, External
Climate Class	SN
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	45
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	Yes